




Sample name: M20966-Ligandrol (LGD-4033, VK-5211)-01 **Injection date:**
Instrument: A010 **Location:**
Inj. volume: 3 **Type:**
Acq. method: HPLC **Sample amount:**
Processing method:
Manually modified:



Signal: DAD1G, Sig=280

RT [min]	Type	Width [min]	Area	Height	Area%
7.909	MM m	0.28	1.59	0.30	0.01
8.217	MM m	0.23	2.32	0.36	0.02
8.716	MM m	0.28	2.21	0.43	0.02
8.874	MM m	0.12	1.18	0.28	0.01
9.029	MM m	0.20	4.59	0.86	0.04
9.333	MM m	0.34	2.92	0.43	0.02
10.161	MM m	0.77	12299.69	1687.13	98.89
10.565	MM m	0.46	121.26	21.25	0.97
12.576	MM m	0.21	2.50	0.56	0.02