

# Shimadzu LCMS-029 Report

## Sample Information

LCMS

Acquired by : System Administrator

Sample Name : M2104-BYL719-06-1

Method File : B 5%-95%-P-new - 2.6min.lcm

Injection Volume : 4

Vial : 21


Column : Waters SunFire C18 50\*4.6mm 5um 2.000ml/min 2.6min

Pump A : 0.1% FA and 10% ACN in H2O


Pump B : 0.1% FA and 10% H2O in ACN

## Chromatogram


mAU



mAU



段号1



## Integration Results for DAD

### PDA Ch1 254nm

| Peak# | Ret. Time | Height | Area    | Area%   |
|-------|-----------|--------|---------|---------|
| 1     | 1.553     | 753188 | 1511416 | 100.000 |
| 总计    |           | 753188 | 1511416 | 100.000 |


### PDA Ch2 214nm

| Peak# | Ret. Time | Height  | Area    | Area%   |
|-------|-----------|---------|---------|---------|
| 1     | 1.553     | 1535850 | 3948385 | 100.000 |


| Peak# | Ret. Time | Height  | Area    | Area%   |
|-------|-----------|---------|---------|---------|
| 总计    |           | 1535850 | 3948385 | 100.000 |

# Mass Spectrum


Rt=1.182min 正




Rt=1.345min 负




Rt=1.362min 正




Rt=1.502min 负




Rt=1.519min 正




Rt=2.065min 负




Rt=2.104min 正



Rt=2.245min 负



Rt=2.402min 负



Rt=2.397min 正

