



Agilent Technologies

Sample name: M21084-GA-017-02

Instrument: A010

Inj. volume: 1.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1A, Sig=220.0, 4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
5.330		4.82	5.6	0.120
6.288		59.11	40.0	1.471
7.221		3945.30	747.5	98.177
8.471		9.32	6.4	0.232