



Agilent Technologies

Sample name: M2112-LY2835219 mesylate-05

Injection date:

Instrument: A010

Location:

Inj. volume: 3


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1C,Sig=254,4 Ref=off

RT [min]	Peak Resolution USP	Width [min]	Area	Height	Area%
6.025		0.11	5.1	1.6	0.04
6.198		0.52	13958.8	3458.3	99.95
6.650		0.13	2.5	0.7	0.02