




LCMS Report

Sample Name : M2127-Hesperadin-01
Method File : B 5%-95%-P.lcm

Method info. : Waters SunFire C18 50*4.6mm 5um 2.000ml/min 2.6min Column Temperature: 40 °C
Gradient: 5% B hold for 0.2min, increase to 95 % B within 1.40 min, hold at 95 % B for 0.9 min, then back to 5% B within 0.01 min

Pump A : 0.03% TFA in H2O
Pump B : 0.03% TFA in ACN

Chromatogram



Integration Results for DAD

PDA Ch1 254nm

Peak#	Ret. Time	Height	Area	Area%
1	1.229	415685	788681	100.000
Total		415685	788681	100.000

PDA Ch2 214nm

Peak#	Ret. Time	Height	Area	Area%
1	1.229	1180215	2257330	100.000
Total		1180215	2257330	100.000

Mass Spectrum

Rt=1.229min Positive

