



Agilent Technologies

Sample name: M21575-Ervogastat-01

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1A,Sig=220.0,4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
6.039		1.9	0.4	0.023
6.141		0.4	0.2	0.004
6.519		8476.4	2156.9	99.923
6.831		4.2	1.0	0.049