




Sample name: M21648-PDK4-IN-1 hydrochloride-01 Injection date:
Instrument: 1 Location:
Inj. volume: HPLC Type:
Acq. method: Sample amount:
Processing method:
Manually modified:



Signal: DAD1C,Sig=254,4 Ref=off

RT [min]	Peak Resolution USP	Width [min]	Area	Height	Area%
6.150		0.05	0.6	0.2	0.01
6.541		0.06	6933.2	1869.4	99.68
6.859		0.05	3.3	1.0	0.05
7.294		0.05	4.5	1.3	0.06
7.514		0.08	1.5	0.3	0.02
7.981		0.05	1.5	0.5	0.02
8.073		0.04	0.4	0.2	0.01
8.331		0.06	2.5	0.7	0.04
8.445		0.06	2.0	0.5	0.03
8.768		0.07	4.0	0.9	0.06
9.074		0.06	1.8	0.5	0.03
	Sum		6955.37		