



Agilent Technologies

Sample name: M21680-TJ-M2010-5-01

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.677	BB	0.58	11010.98	2873.73	99.95
8.899	MB m	0.15	5.93	1.49	0.05