



8.326 8.306 7.788 7.771 7.767 7.750 7.731 7.672 7.667 7.651 7.646 7.557 7.538 7.528 7.522 7.503 7.486 7.482 7.426 7.408 7.388 7.337 7.319 7.300 7.260 7.219 7.196 7.133 7.111 6.822 6.801 4.600 4.586 4.570 4.570 3.415 3.347 3.299 3.294 2.491 2.463 2.245 2.230 2.204 2.185 2.171 1.666 1.255 1.241 1.191 1.175 -0.002



F1: 399.751	F2: 100.526	SW1: 8013	OF1: 3200.8	PTS1d: 16384
EX: s2pul	PW: 5.9 usec	PD: 1.0 sec	NA: 127	LB: 0.0
Sample Name: M21708-Erastin2-01				