




Sample name: M2180-SGI-7079-02  
Instrument: A012  
Inj. volume: 1.5  
Acq. method: HPLC  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1B,Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Area%
5.966	11.08	3.80	0.16
6.221	6813.99	2167.16	99.22
6.408	12.91	3.17	0.19
6.634	3.86	1.39	0.06
7.865	25.42	8.40	0.37