




Sample name: M2206-TH-302-05
Instrument: 4.5
Inj. volume: HPLC
Acq. method:
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1B, Sig=210,4 Ref=off

RT [min]	Area	Height	Area%
6.255	25.85	8.26	0.57
6.333	8.93	3.32	0.20
6.624	40.38	11.08	0.88
6.751	6.40	1.97	0.14
7.009	4488.80	1170.09	98.22
Sum	4570.37		