



Agilent Technologies

Sample name: M22523-Recilisib sodium-02

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1A, Sig=220.0, 4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
7.389		33.3	9.9	0.222
7.929		223.3	53.7	1.489
8.208		20.0	6.2	0.134
8.583		14719.3	2806.1	98.156