



Agilent Technologies

Sample name: M2293-CX-5461-04

Instrument: A010

Inj. volume: 1.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1E,Sig=240.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.123	BB	0.51	5962.48	777.27	99.50
6.564	BB	0.22	11.97	2.64	0.20
7.540	BM m	0.13	12.73	3.06	0.21
7.607	MB m	0.17	5.09	1.02	0.08