



Sample name: M2303-Capivasertib (AZD5363)-04

Injection date:

Instrument: A012

Location:

Inj. volume: 2.5


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1A,Sig=220.0,4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
4.923		6.7	2.3	0.072
5.146		16.6	3.8	0.179
5.412		10.4	5.7	0.113
5.491		9180.8	2515.2	99.255
5.927		15.6	3.6	0.169
6.063		3.3	0.9	0.036
6.303		2.0	0.7	0.022
6.577		14.2	4.6	0.154