



Agilent Technologies

Sample name: M2376-Aloxistatin-02

Instrument: A012

Inj. volume: 8

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B,Sig=210,4 Ref=off

RT [min]	Area	Height	Area%
4.795	3.08	1.27	0.03
4.939	3.60	1.21	0.04
5.156	5.90	1.45	0.06
7.245	29.60	8.87	0.31
8.059	10.67	3.02	0.11
8.590	9509.73	1544.40	98.53
8.922	4.67	1.15	0.05
9.116	7.56	1.12	0.08
9.216	3.90	1.08	0.04
9.857	72.60	13.97	0.75