



Agilent Technologies

Sample name:

M25404-Nerandomilast (BI 1015550)-02

Instrument:

A011

Inj. volume:

4

Acq. method:

HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1E, Sig=240,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
5.852	BB	0.60	15570.68	3262.25	99.99
7.465	MM m	0.09	1.26	0.42	0.01
7.604	MM m	0.10	1.04	0.34	0.01