



Agilent Technologies

Sample name: M25459-Zunsemetinib-01-1

Injection date:

Instrument: A012

Location:

Inj. volume: 2


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1A, Sig=220.0,4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
6.556		9.1	3.0	0.088
7.363		7.5	2.2	0.072
7.600		10323.5	2555.8	99.421
8.116		36.8	9.4	0.354
9.959		6.8	1.4	0.065