


SAMPLE INFORMATION

Sample Name:	M25556-SHA 68-04	Acquired By:	System
Sample Type:	Unknown	Sample Set Name:	M25556
Vial:	85	Acq. Method Set:	CD_05_95_15min_40D
Injection #:	1	Processing Method:	LC PQ
Injection Volume:	6.00 ul	Channel Name:	220.0nm
Run Time:	15.0 Minutes	Proc. Chnl. Descr.:	W2996 PDA 220.0 nm (PDA)



	RT	Area	% Area	Height
1	3.647	28453	0.30	5595
2	5.001	15102	0.16	3377
3	6.857	11406	0.12	3101
4	7.889	13566	0.14	3964
5	8.637	9534945	99.03	1917561
6	9.263	24896	0.26	5928