



Agilent Technologies

Sample name: M2558-DMOG-06

Instrument: A010

Inj. volume: 4

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1A, Sig=220.0, 4.0 Ref=off

RT [min]	Peak	Area	Height	Area%
3.420		15856.8	1387.8	99.534
4.338		66.2	12.4	0.416
7.275		2.5	0.7	0.016
9.384		5.5	1.2	0.035