



Sample name: M2577-StemRegenin 1-02

Injection date:

Instrument: 5

Location:

Inj. volume: HPLC


Type:

Acq. method:

Sample amount:

Processing method:

Manually modified:



Signal: DAD1A,Sig=220,4 Ref=off

RT [min]	Peak Resolution USP	Width [min]	Area	Height	Area%
7.466		0.22	32.7	7.9	0.185
8.615		0.06	2.1	0.6	0.012
8.839		0.07	5.1	1.3	0.029
9.413		0.60	17577.1	2988.2	99.521
10.631		0.21	31.6	6.4	0.179
16.064		0.08	13.1	2.7	0.074
Sum			17661.71		