



Agilent Technologies

Sample name: M2731-GSK-J4 hydrochloride-06

Injection date:

Instrument: A010

Location:

Inj. volume: 2


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1A,Sig=220.0,4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
5.709		3.9	1.3	0.033
5.808		2.3	0.9	0.019
5.957		14.0	4.6	0.118
6.055		4.5	1.6	0.038
6.112		6.8	2.4	0.057
6.483		2.0	0.5	0.017
7.429		5.4	1.3	0.046
7.601		4.0	1.4	0.033
7.853		11716.4	2707.8	98.955
7.981		29.7	9.7	0.250
8.098		10.5	2.2	0.088
8.387		16.3	3.6	0.138
8.756		3.3	0.9	0.028
9.033		4.8	1.1	0.040
9.370		16.4	3.3	0.138