




Sample name: M2784-JNJ-7777120-02
Instrument: 2
Inj. volume: HPLC
Acq. method:
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1A,Sig=220,4 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
5.650		4.2	1.2	0.032
5.905		4.8	1.6	0.037
6.055		12810.3	2920.7	98.818
6.273		4.3	1.7	0.033
6.508		1.7	0.6	0.013
6.723		3.8	1.0	0.029
6.865		0.9	0.4	0.007
6.953		3.8	1.2	0.029
7.928		8.4	2.1	0.065
8.171		59.1	17.0	0.456
9.449		7.8	2.1	0.060
9.685		14.1	3.3	0.109
10.173		5.1	1.2	0.039
10.450		35.5	7.4	0.274
	Sum	12963.50		