


Data file: C:\CHEM32\1\DATA\240813\M28234-03Final.D
Sample name: M28234- Hepcidin-25 (human)-03
Description: Mobile Phase: Buffer A: 0.1% TFA in H₂O
Buffer B: 0.1% TFA in ACN
Gradient: 24-34%B in 20min, Flow rate: 1.0ml/min 50C
Column: Unitary C18 5um 100A 4.6*150mm
Instrument: Agilent HPLC 1100

Instrument: 1100
Injection date: 2024/8/13 11:54:58
Last changed: 2024/8/13 13:02:40

Location: 75
Injection: 1 of 1
Injection volume: 20.000
Acq.operator: SYSTEM



Signal: DAD1 A, Sig=220, 4 Ref=360, 100

RT[min]	Type	Width[min]	Area	Height	Area%	Name
9.234	MF	0.5477	8080.7476	245.8957	98.1593	
10.145	FM	0.4755	151.5300	5.3108	1.8407	
Sum			8232.2775			