



Sample name: M28266-KL-11743-01

Instrument: A010

Inj. volume: 1.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B, Sig=210,4 Ref=off

RT [min]	Area	Height	Area%
6.772	121.67	21.36	1.00
7.087	11946.06	1701.37	98.10
7.395	91.03	28.33	0.75
8.004	5.99	2.00	0.05
8.147	6.00	1.31	0.05
8.642	6.21	1.79	0.05