


SAMPLE INFORMATION

Sample Name:	M28606-CPUY192018-02	Acquired By:	System
Sample Type:	Unknown	Sample Set Name:	M28606
Vial:	82	Acq. Method Set:	CD_05_95_15min_40D
Injection #:	1	Processing Method:	LC PQ
Injection Volume:	4.00 ul	Channel Name:	235.0nm
Run Time:	15.0 Minutes	Proc. Chnl. Descr.:	W2996 PDA 235.0 nm (PDA)



	RT	Area	% Area	Height
1	6.843	7627	0.06	2164
2	7.210	56116	0.44	13478
3	7.965	12525049	99.05	2562032
4	8.427	55271	0.44	13338
5	9.012	1352	0.01	335