



Agilent Technologies

Sample name: M28813-Ro 25-6981-02

Instrument: A010

Inj. volume: 4

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B, Sig=210.0, 4.0 Ref=off

RT [min]	Area	Height	Area%
6.189	48.36	16.31	0.57
6.936	8424.41	2004.37	99.27
7.262	3.59	0.99	0.04
7.481	4.11	1.49	0.05
8.010	5.73	1.83	0.07