


Data file: C:\CHEM32\1\DATA\231201\M29001Final.D
 Sample name: M29001-Bremelanotide-02
 Description: Mobile Phase: Buffer A: 0.1% TFA in H2O
 Buffer B: 0.1% TFA in ACN
 Gradient: 27-37% in 20min, Flow rate: 1.0ml/min
 Column: Unitary C18 5um 100A 4.6*150mm
 Instrument: Agilent HPLC 1100

Instrument: 1100
 Injection date: 2023/12/1 15:53:18
 Last changed: 2023/12/1 16:29:34

Location: 31
 Injection: 1 of 1
 Injection volume: 10.000
 Acq.operator:



Signal: DAD1 A, Sig=220,4 Ref=360,100

RT[min]	Type	Width[min]	Area	Height	Area%	Name
6.825	MM	0.1510	3.4194	0.3775	0.0252	
7.124	MM	0.1613	3.1594	0.3264	0.0233	
7.664	BV F	0.1554	11.3722	1.1497	0.0837	
7.845	VB	0.1426	10.1518	1.1544	0.0747	
8.453	VV E	0.1853	71.8885	6.1064	0.5292	
8.872	VBAR	0.2266	13485.3174	967.0767	99.2640	
Sum			13585.3088			