




Sample name: M2960-Probenecid-04  
Instrument: A010  
Inj. volume: 3  
Acq. method: HPLC  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1C,Sig=254.0,4.0 Ref=off

RT [min]	Peak	Width [min]	Area	Height	Area%
7.009		0.22	33.7	10.4	0.35
7.844		0.11	0.4	0.1	0.00
8.946		0.73	9609.6	2007.5	99.64
13.162		0.17	0.9	0.2	0.01