




Agilent Technologies

Sample name: M2964-Propranolol hydrochloride-03 Injection date:
Instrument: A010 Location:
Inj. volume: 2.5 Type:
Acq. method: HPLC Sample amount:
Processing method:
Manually modified:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
4.982	MM m	0.10	0.50	0.18	0.00
5.084	MM m	0.10	2.39	0.88	0.01
6.264	MB m	0.12	5.46	1.69	0.03
6.381	MM m	0.09	1.24	0.48	0.01
6.587	BB	0.09	3.35	1.30	0.02
6.752	BB	0.45	17632.33	2949.09	99.82
9.014	MM m	0.15	14.35	3.23	0.08
9.074	MB m	0.08	4.29	1.56	0.02