


Data file: D:\CHEM32\1\DATA\24110201Final.D
 Sample name: M2976-RGD-04
 Description: Mobile phase: Buffer A: 0.1%TFA in H2O
 Buffer B: 0.1%TFA in ACN
 Gradient: 0-0% B in 20min , Flow rate: 1.0ml/min
 Column: Unitary C18, 4.6*150mm, 5um 100A
 Instrument: Agilent HPLC 1100

Instrument: 1100

Location: 61
 Injection: 1 of 1
 Injection volume: 10.000
 Acq.operator: SYSTEM



Signal: DAD1 A, Sig=220,4 Ref=360,100

RT[min]	Type	Width[min]	Area	Height	Area%	Name
2.975	BB	0.2392	5.1616	0.3602	0.1417	
3.569	BB	0.1384	2.1868	0.2592	0.0600	
4.018	BV E	0.1633	33.9811	3.4394	0.9327	
4.391	VV R	0.3124	3582.4353	162.9774	98.3241	
6.200	VB E	0.2954	19.7316	1.0306	0.5416	
Sum			3643.4963			