




Sample name: M29821-HS-276-01
Instrument: A010
Inj. volume: 5
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1B, Sig=210.0, 4.0 Ref=off

RT [min]	Area	Height	Area%
5.344	6.44	1.96	0.10
5.604	2.46	0.85	0.04
5.742	6354.63	1552.19	98.52
7.311	23.44	7.57	0.36
7.384	62.90	20.34	0.98