



Agilent Technologies

Sample name: M3001-SC144-01

Instrument: A012

Inj. volume: 1.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1A,Sig=220,4 Ref=off

RT [min]	Peak	Area	Height	Area%
5.426		6932.6	2159.8	99.861
6.229		3.6	1.1	0.051
6.357		1.3	0.5	0.019
7.284		1.5	0.6	0.021
8.138		3.3	1.0	0.048