




Agilent Technologies

Sample name: M30583-8-OH-DPAT hydrobromide-01 Injection date:
Instrument: A010 Location:
Inj. volume: 3 Type:
Acq. method: HPLC Sample amount:
Processing method:
Manually modified:



Signal: DAD1B,Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Area%
5.597	2.81	0.92	0.05
6.410	5467.96	1405.03	99.56
6.657	2.79	0.80	0.05
6.825	2.99	1.07	0.05
6.944	12.75	4.24	0.23
7.085	2.86	0.95	0.05