

Single Injection Report



Agilent Technologies

Data file: 02.dx

Sequence Name: 20240315-6

Sample name: M31311-Tesamorelin-03

Instrument: SB-163

Inj. volume: 5.000

Acq. method: YGWZ09.amx

Processing method: YGWZ.pmx

Manually modified: Manual Integration

Project Name: CRO


Operator: lanxinling

Injection date: 2024-03-15 20:34:23+08:00

Location: 37

Type: Sample

Sample amount: 0.00



Signal: DAD1A,Sig=220,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%	Name
8.512	MM m	0.10	5.48	0.71	0.10	
8.700	MM m	0.14	3.60	0.32	0.07	
9.086	MM m	0.07	1.58	0.33	0.03	
9.205	MM m	0.12	5.44	0.59	0.10	
10.466	MM m	0.09	2.98	0.47	0.06	
10.608	MM m	0.12	3.64	0.40	0.07	
10.895	MM m	0.04	1.63	0.66	0.03	
11.135	MM m	0.24	5279.18	346.11	98.98	
11.575	MM m	0.14	25.03	3.00	0.47	
11.755	MM m	0.08	5.11	1.02	0.10	
Sum			5333.66			