

Data file: C:\CHEM32\1\DATA\240125\24012505Final.D
Sample name: M31313-Pegcetacoplan
Description: Mobile Phase: Buffer A: 0.1% TFA in H₂O
Buffer B: 0.1% TFA in ACN
Gradient: 0-100%B in 30min, Flow rate: 1.0ml/min
Column: UniPS 10-300 10um 100A 4.6*250mm
Instrument: Agilent HPLC 1100

Instrument: 1100

Injection date: 2024/1/25 14:06:58


Last changed: 2024/1/25 14:25:14

Location: 91

Injection: 1 of 1

Injection volume: 10.000

Acq.operator: System



Signal: DAD1 B, Sig=280,4 Ref=off

RT[min]	Type	Width[min]	Area	Height	Area%	Name
13.879	BV	0.2285	5.1202	0.3462	0.4430	
15.191	BB	0.1367	2.7234	0.3284	0.2356	
16.874	BB	0.1898	1141.1255	99.4878	98.7299	
22.606	BB	0.3017	6.8365	0.3356	0.5915	
Sum			1155.8056			