


Sample Name: M3186-Isoprenaline hydrochloride-04-1

=====
Acq. Operator : Z
Acq. Instrument : instr1

Location : Vial 51

Inj Volume : 5.000 µl

Acq. Method : D:\CHEM32\1\METHODS\Z2.M



=====
Area Percent Report
=====

Sorted By : Signal
Multiplier: : 1.0000
Dilution: : 1.0000
Use Multiplier & Dilution Factor with ISTDs

Signal 1: VWD1 A, Wavelength=220 nm

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	4.154	BV	0.1875	86.84724	6.43344	0.3725
2	4.726	VB	0.1348	2.31678e4	2604.77734	99.3760
3	6.759	BB	0.0710	58.62274	12.21678	0.2515

Totals : 2.33133e4 2623.42757

=====
*** End of Report ***
=====