

Sample Name : M3192-TRAM-34-01

Method Info : Mobile Phase : A:water(0.01%TFA) B:ACN(0.01%TFA)

Gradient:5% increase to 95%B within 1.3min,95%B for 1.7min


Flow Rate :2ml/min

Column :Sunfire, 50*4.6mm,3.5um


Column Temperature:50C

Detection:UV(214,4nm) and MS(ESI,Pos mode,110 to 1000 amu)


*ADC1 A, ELSD




*DAD1 A, Sig=2144 Ref=off



*DAD1 B, Sig=2544 Ref=off



MSD1 TIC, MS File ES-AP I, Pos, Scan, Frag: 70



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
2.18	0.02	96.33	95.62	100.00

Integration Results for DAD1 A, Sig=2144 Ref=off

RetTim	Width	Area	Height	Area%
2.13	0.02	16.83	12.61	0.45
2.18	0.03	3716.40	2074.76	99.55

Integration Results for DAD1 B, Sig=2544 Ref=off

RetTim	Width	Area	Height	Area%
2.18	0.02	79.60	58.51	100.00

Ret. Time: 2.18 <<<< POSITIVE SPECTRA >>>>

