



Agilent Technologies

Sample name: M3444-A-769662-03

Instrument: A012

Inj. volume: 3

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1C,Sig=254,4 Ref=off

RT [min]	Peak	Width [min]	Area	Height	Area%
6.212		0.15	5.3	1.7	0.04
6.384		0.11	0.5	0.2	0.00
6.543		0.11	2.1	0.7	0.02
6.652		0.16	1.6	0.4	0.01
6.937		0.26	30.7	8.6	0.24
7.313		0.16	2.9	0.7	0.02
7.520		0.13	0.8	0.2	0.01
7.786		0.84	12545.2	3126.5	99.54
9.440		0.16	5.1	1.3	0.04
9.760		0.36	8.8	1.3	0.07