



Agilent Technologies

Sample name: M3508-Triamcinolone Acetonide-02

Injection date:

Instrument: A010

Location:

Inj. volume: 2


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1F,Sig=240.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.136	MM m	0.13	2.40	0.56	0.04
6.472	MM m	0.07	0.63	0.25	0.01
7.314	BB	0.16	8.89	2.62	0.15
7.529	MM m	0.06	1.16	0.56	0.02
7.670	MB m	0.44	6109.12	1792.64	99.73
8.022	MM m	0.10	1.54	0.54	0.03
8.918	MM m	0.12	2.14	0.60	0.03