



Agilent Technologies

Sample name: M3666-Sofosbuvir (PSI-7977)-02

Injection date:

Instrument: A012

Location:

Inj. volume: 2


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1E, Sig=260.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
7.387	BM m	0.49	3828.56	1185.14	99.96
7.885	MM m	0.21	1.63	0.31	0.04