



Sample name: M3677-N6022-02

Instrument: A012

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
5.571	BM m	0.26	9294.79	2542.78	99.13
5.745	MB m	0.12	5.41	1.59	0.06
5.971	BM m	0.10	4.93	1.57	0.05
6.023	MB m	0.12	30.17	8.69	0.32
6.144	BB	0.15	25.45	8.48	0.27
6.392	MB m	0.11	13.56	4.92	0.14
6.791	MM m	0.09	2.19	0.78	0.02