



Agilent Technologies

Sample name: M3699-Letrozole-02

Instrument: A012

Inj. volume: 1

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.819	MM m	0.08	1.22	0.43	0.01
7.693	BB	0.60	15155.80	3073.54	99.96
8.152	BB	0.13	5.38	1.68	0.04