



Agilent Technologies

Sample name: M3864-Remdesivir-02

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1E, Sig=240.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
5.901	MM m	0.15	1.14	0.30	0.03
6.218	MM m	0.15	1.06	0.23	0.02
6.681	MM m	0.13	1.15	0.33	0.03
6.975	MM m	0.08	1.15	0.46	0.03
7.083	MB m	0.11	4.84	1.64	0.11
7.602	BM m	0.31	4417.09	1318.63	99.66
7.739	MM m	0.06	3.31	1.32	0.07
7.838	MM m	0.07	1.30	0.37	0.03
8.595	MM m	0.13	1.29	0.32	0.03