

Sample Name : Celastrol-04




Method Info : Column:SunFire C18 (4.6x 50 mm,3.5 um)

Mobile phase: H2O(0.01%TFA) (A) / ACN(0.01%TFA) (B)

Elution program: Gradient from 10 to 95% of B in 1.4 min at 2ml/min

Temperature: 50 °C

Detection: UV (214 , 4 nm) and MS (ESI, Pos mode ,110 to 1000 amu)



Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
2.27	0.04	378.98	155.40	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
2.27	0.04	280.19	116.52	100.00

Ret. Time: 2.14

<<<< POSITIVE SPECTRA >>>>



Ret. Time: 2.27

