




Sample name: M3970-Scutellarin-03  
Instrument: A010  
Inj. volume: 5  
Acq. method: HPLC  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1B, Sig=210.0, 4.0 Ref=off

RT [min]	Area	Height	Area%
5.135	4.39	1.44	0.11
5.510	4035.45	824.66	98.10
5.913	60.58	11.29	1.47
6.483	7.27	1.25	0.18
6.925	5.95	1.42	0.14