



Agilent Technologies

Sample name: M3973-Tazarotene-01

Instrument: 2

Inj. volume: HPLC

Acq. method:

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1C,Sig=254,4 Ref=off

RT [min]	Width [min]	Area	Height	Area%
10.563	0.28	2.8	0.5	0.05
11.076	0.17	0.9	0.2	0.02
11.482	0.21	0.9	0.2	0.01
14.380	1.24	6059.0	657.5	99.90
15.088	0.18	1.7	0.4	0.03