



Agilent Technologies

Sample name: M3975-Lorlatinib (PF-06463922)-02

Injection date:

Instrument: A012

Location:

Inj. volume: 4


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
5.580	MM m	0.09	1.84	0.67	0.02
5.720	MM m	0.06	0.18	0.10	0.00
5.857	BB	0.41	11968.35	3074.83	99.52
6.273	BB	0.12	7.93	2.78	0.07
6.739	MB m	0.11	8.50	2.81	0.07
6.911	MM m	0.09	1.40	0.51	0.01
7.141	BB	0.12	9.51	3.20	0.08
7.487	MM m	0.10	2.25	0.80	0.02
7.649	MB m	0.14	23.28	7.42	0.19
7.865	MB m	0.09	2.71	0.96	0.02