




Sample name: M3994-Hederagenin-03
Instrument: A010
Inj. volume: 10
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1H,Sig=200,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
11.765	MB m	0.33	1513.36	252.75	98.87
11.921	BB	0.07	8.78	3.60	0.57
11.995	BB	0.07	8.44	3.14	0.55