



Agilent Technologies

Sample name: M3998-Rhein-03

Instrument: A012

Inj. volume: 10

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
7.783	MM m	0.07	0.64	0.21	0.03
8.221	MB m	0.11	3.37	1.05	0.14
8.614	MM m	0.12	3.07	0.89	0.13
8.851	BM m	0.61	2385.00	482.67	99.38
9.861	MB m	0.17	7.78	1.76	0.32