




Sample name: M40519-AOH1996-02  
Instrument: A010  
Inj. volume: 1  
Acq. method: HPLC  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1A, Sig=220.0, 4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
7.490		5.2	1.6	0.032
8.071		1.0	0.3	0.006
8.526		184.1	50.7	1.140
10.588		15856.7	2436.1	98.176
10.883		72.1	16.8	0.446
10.998		0.7	0.3	0.004
11.259		13.2	3.5	0.082
11.358		3.5	1.1	0.022
11.588		2.0	0.6	0.013
11.804		6.2	1.6	0.038
11.944		6.7	1.7	0.041